MicroLife WS 80

- 1) ONSET Button
2) User Button
3) Down Button
4) Up Button
5) Battery Compartment
6) KG / LB / ST Switch

- Display
7) Low Battery Indication
8) Zero Setting
9) Weight Reading
10) Gender
11) Age
12) Height
13) Body Fat Reading
14) Body Water Reading
15) Muscle Mass Reading
16) Bone Mass Reading
17) Body Mass Index (BMI) Reading
18) Basal Metabolic Rate (BMR) Reading
19) Error

Read the instructions carefully before using this device.

Dear Customer
This diagnostic scale allows you to evaluate your body composition. The knowledge of your body composition allows you to combine a correct nutrition and physical activity as essential conditions to improve the quality of life, promote health, and prevent pathologies.

1. Important Facts about Body Composition

It is important to know the body composition of an individual in order to be able to estimate his/her health condition. It is not enough to know one's weight but it is necessary to identify its quality, which means dividing the human body into its components...

- Lean body mass: bones, muscles, viscera, water, organs, blood
Plyometric mass: adipose cells

Body Fat
The total body fat mass is calculated by totalling the essential fat (or primary fat) and the spare fat (or storage fat). The essential fat is necessary for the metabolism of some spare fat...

Body Water
Body water % is an important indication of the overall well-being of a human being. It is the single most important component of body weight. Research shows that 57 % body water is optimum...

Muscle Mass
Muscle mass is important in determining a healthy body composition. A person with a higher % of muscle mass finds it easier to move, but needs more energy to do it. Exercise is very important in maintaining a healthy body composition...

Bone Mass
Bone mass increases rapidly in childhood and reaches its maximum between 30-40 years. It decreases slightly afterwards with increasing age. Factors like age, gender, weight, and height have a small influence on the bone mass...

Body Mass Index (BMI)
The BMI is a simple index of weight-for-height that is commonly used to classify underweight, overweight and obesity in adults. It is calculated by dividing the body weight by the square of the height in m of a person.

Basal Metabolic Rate (BMR)
The BMR is kcal is a measure of energy required by the body when in a state of complete rest to maintain its basic functions. About 70 % of a human's total energy expenditure is due to the basal metabolic rate...

2. Suggestions for Use
Placing the scale on a flat, hard surface. Soft, inaccurate floors (e.g. rugs, carpets, linoleum) can cause the scale to give uneven readings...

3. Using the Device for the First Time
How to insert batteries. The scale is supplied with two AA batteries. The battery compartment is located on the bottom of the scale...

4. Entering Personal Data
This scale can memorize the data of 10 different persons. Turn the scale on by pressing the ONSET button. Turn the scale on by pressing the ONSET button...

5. Calculation of Body Composition
Turn the scale on by pressing the ONSET button. Press the 'C' symbol. Select your memory location and press USER (2) to confirm. The data entered as well as the last measured values will be viewed on the display...

6. Important Safety Instructions
This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application...

7. Garantie
This device is covered by a 2 year guarantee period. The date of purchase. During this guarantee period, at our discretion, MicroLife will repair or replace the defective product free of charge...

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Muscle mass
Muscle mass is important for determining a healthy body composition. A person with a higher % of muscle mass finds it easier to move, but needs more energy to do it. Exercise is very important in maintaining a healthy body composition...

Bone mass
Bone mass increases rapidly in childhood and reaches its maximum between 30-40 years. It decreases slightly afterwards with increasing age. Factors like age, gender, weight, and height have a small influence on the bone mass...

Body mass index (BMI)
The BMI is a simple index of weight-for-height that is commonly used to classify underweight, overweight and obesity in adults. It is calculated by dividing the body weight by the square of the height in m of a person.

Basal metabolic rate (BMR)
The BMR is kcal is a measure of energy required by the body when in a state of complete rest to maintain its basic functions. About 70 % of a human's total energy expenditure is due to the basal metabolic rate...

2. Suggestions for use
Placing the scale on a flat, hard surface. Soft, inaccurate floors (e.g. rugs, carpets, linoleum) can cause the scale to give uneven readings...

3. Using the device for the first time
How to insert batteries. The scale is supplied with two AA batteries. The battery compartment is located on the bottom of the scale...

4. Entering personal data
This scale can memorize the data of 10 different persons. Turn the scale on by pressing the ONSET button. Turn the scale on by pressing the ONSET button...

5. Calculation of body composition
Turn the scale on by pressing the ONSET button. Press the 'C' symbol. Select your memory location and press USER (2) to confirm. The data entered as well as the last measured values will be viewed on the display...

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