

Demo APP for MicroLifeDeviceSDK (WatchBP Office) - iOS


Table of Contents

- Chapter 1 Development Environment**
- Chapter 2 Entry Point and Bluetooth LE Protocol
 (BLEManager)**
- Chapter 3 APIs of BLEManager**
- Chapter 4 WatchBP Office APIs**
- Chapter 5 User Interface & Functionality of WatchBP Office**
- Chapter 6 Demo App**

Chapter 1 Development Environment

This user manual serves as a quick introduction to MicroLifeDeviceSDK / APIs and shows how to integrate into an iOS Demo App.

- 1.1. First of all, add MicroLifeDeviceSDK.framework into a development project. The minimum supported version is iOS 10. The compatible version of Xcode IDE is 10.2.1.
- 1.2. Under TARGETS / General / Linked Frameworks and Libraries, add MicroLifeDeviceSDK.framework.



- 1.3. Under TARGETS / Build Settings / Enable, set Bitcode to NO.
- 1.4. Under TARGETS / Build Settings / Other Linker Flags, use the linker flag “-all_load”.
- 1.5. Import header file as bellows.

```
#import <MicroLifeDeviceSDK/MicroLifeDeviceSDK.h>
```

Chapter 2 Entry Point and Bluetooth LE Protocol (BLEManager)

The BLEManager is responsible for managing the Bluetooth communication. Refer to the file

“MicroLifeDeviceSDK/MicroLifeDeviceSDK.h” for suitable device models.

They are named like “XXXBLEManager”.



- 2.1. Initiate a BLEManager and register the delegation of BLEManager. The DataResponseDelegate is dedicated for response messages of BLEManager.
- 2.2. Once the BLEManager is initiated, the BLEManager checks the status of Bluetooth by
BLEManagerCellPhoneBluetoothDidUpdateState.
- 2.3. While turning on Bluetooth, BLEManager runs the scan/discovery procedure automatically by
BLEManagerDidDiscoverBluetoothDeviceMacAddress to get devices' information in the vicinity.
- 2.4. Call the bindingDevice:(NSData *)macAddress to bind an activated device. The status of binding can be found by the following section 2.5. The BLEManager will automatically connect to a bound device while it is activated.
- 2.5. Connection status :
 - 2.5.1. Connection / binding :

- (void) * BLEManagerDidConnectDevice;

2.5.2. Fail of connection :

- (void) * BLEManagerDidFailToConnectDevice;

2.5.3. Disconnection :

- (void) * BLEManagerDidDisconnectDevice;

Chapter 3 APIs of BLEManager

By utilizing the APIs of BLEManager, can transfer data between device and App via Bluetooth. The APIs includes two parts: one is for fundamental Bluetooth and another is dedicated for communication with a designated device / BPM (Blood Pressure Monitor) with specific commands and protocol.

3.1. Device searching :

3.1.1.

	<code>+(instancetype)shareInstanceWhithAuthorizationkey:(NSString *)key Target * Names:(NSArray *)target * Names</code>
Definition	Searching for device with default / customized name
Parameter	Key: Authorization code target * Names: device name which can be included single or multiple
	<pre>self.aBPMBLEManager = [BPMBLEManager shareInstanceWhithAuthorizationkey:\$DKkey TargetBPMNames:@[@"A6 BT",@"A6 BASIS PLUS BT",@"A7 TCUCH BT",@"B3 BT",@"B6 Connect",@"A6",@"Progress"]]; self.aBPMBLEManager.dataResponseDelegate = self;</pre>

3.1.2.

	<code>-(void) *</code> <code>BLEManagerCellPhoneBluetoothDidUpdateState:(MicroLifeBLEState)state;</code>
Definition	Response for the shareInstanceWhithAuthorizationkey (). This is to manage Bluetooth status of cellphone.
Parameter	status: Bluetooth status of cellphone

	<pre> * @enum MicroLifeBLEState * * @discussion Represents the current state of a BLE. * * @constant MicroLifeBLEStateUnknown=CBManagerStateUnknown State unknown, update imminent. * @constant MicroLifeBLEStateResetting=CBManagerStateResetting The connection with the system service was momentarily lost, update imminent. * @constant MicroLifeBLEStateUnsupported=CBManagerStateUnsupported The platform doesn't support the Bluetooth Low Energy Central/Client role. * @constant MicroLifeBLEStateUnauthorized=CBManagerStateUnauthorized The application is not authorized to use the Bluetooth Low Energy role. * @constant MicroLifeBLEStatePoweredOff=CBManagerStatePoweredOff Bluetooth is currently powered off. * @constant MicroLifeBLEStatePoweredOn=CBManagerStatePoweredOn Bluetooth is currently powered on and available to use. * */ typedef NS_ENUM(NSUInteger, MicroLifeBLEState) { MicroLife3LEStateUnknown = 0, MicroLife3LEStateResetting, MicroLife3LEStateUnsupported, MicroLife3LEStateUnauthorized, MicroLife3LEStatePoweredOff, MicroLife3LEStatePoweredOn, } </pre>
--	---

	<p>- (void) *</p> <p>BLEManagerDidDiscoverBluetoothDeviceMacAddress:(NSData *)macAddress Name:(NSString *)name RSSI:(NSNumber *)RSSI;</p>
Definition	<p>Response for the sharedInstanceWhithAuthorizationkey ().</p> <p>This is to manage information of devices that are discovered in the vicinity.</p>
Parameter	<p>macAddress: MAC address of device</p> <p>name: device name</p> <p>RSSI: RSSI</p>

3.2. Binding / Pairing :

3.2.1.

	- (void)bindingDevice:(NSData *)macAddress
Definition	Binding a specified device by MAC
Parameter	macAddress: MAC address of device

3.2.2.

	- (void) * BLEManagerDidConnectDevice;
Definition	Response for the bindingDevice() with the status Connection

	- (void) * BLEManagerDidFailToConnectDevice
Definition	Response for the bindingDevice() with the status Fail of Connection

	- (void) * BLEManagerDidDisconnectDevice;
Definition	Response for the bindingDevice() with the status Disconnection

3.3. MAC address of binding device :

3.3.1.

	- (NSData *)getBindingDevice;
Definition	Get MAC address of a binding device or nil.

3.4. Unbinding :

3.4.1.

	- (void)unBindingDevice;
Definition	Unbinding device

3.4.2.

	- (void) * BLEManagerDidFailToConnectDevice;
Definition	Response for the unBindingDevice()

3.5. Disconnection :

3.5.1.

	- (void)disconnectDevice;
Definition	Disconnect with device

3.5.2.

	- (void) * BLEManagerDidFailToConnectDevice;
Definition	Response for the disconnectDevice ()

3.6. Searching for devices :

3.6.1.

	- (void)reScan;
Definition	Search for devices repeatedly

3.6.2.

	- (void) * BLEManagerDidDiscoverBluetoothDeviceMacAddress:(NSData *)macAddress Name:(NSString *)name RSSI:(NSNumber *)RSSI;
Definition	This is to get information of devices which are discovered in the vicinity.
Parameter	macAddress: MAC of device name: device name RSSI: RSSI

3.7. Stop searching for devices :

3.7.1.

	- (void)stopScan;
Definition	Stop searching

3.7.2. Null.

Chapter 4 WatchBP Office APIs

4.1. Read all history or current data :

4.1.1. Interface :

	-(void)readAllHistorys;
Definition	Read all history or current data

4.1.2. Delegate :

	- (void)WBOBLEManagerResponseReadAllHistorys:(MicroLifeD Record *)data;
Parameter	data : History data

4.2. Read central BP memory data by index :

4.2.1. Interface :

	- (void)readCBPDataWithIndex:(int)index Dformat:(Dformat)dformat;
Definition	4.2. Read central BP memory data by index
Parameter	index : CBP memory index dformat : Data format NoCBPRaw : No CBP raw data LowCBPRaw : low resolution CBP data FullCBPRaw : full CBP raw data

4.2.2. Delegate :

	- (void)WBOBLEManagerResponseReadCBPData:(MicroLifeCB PdataAndCalCBP *)data IsNoData:(BOOL)isNoData;
Parameter	cRecord : CBP data and CalCBP data isNullData : True or False

4.3. Clear all history data :

4.3.1. Interface :

	- (void)clearAllHistorys;
Definition	Clear all history data

4.3.2. Delegate :

	- (void)WBOBLEManagerResponseClearHistory:(BOOL)isSuccess;
Parameter	isSuccess : True or False

4.4. Disconnect the Bluetooth device :

4.4.1. Interface :

	- (void)disconnect;
Definition	Disconnect the Bluetooth device

4.4.2. Delegate :

	- (void)WBOBLEManagerDidDisconnectDevice;
--	---

4.5. Read user ID and version data

4.5.1. Interface :

	- (void)readUserAndVersionData;
Definition	Read user ID and version data

4.5.2. Delegate :

	- (void)WBOBLEManagerResponseReadUserAndVersionData:(MicroLifeUserInfo *)user VersionData:(MicroLifeDeviceInfo *)verData;
Parameter	user : user ID verData : version data

4.6. Write user ID :

4.6.1. Interface :

	- (void)writeUserID:(NSString *)ID;
Definition	Write user ID to device
Parameter	ID : User ID.

4.6.2. Delegate :

	- (void)WBOBLEManagerResponseWriteUserID:(BOOL)isSuccess;
Parameter	isSuccess : True or False

4.7. Read BPM setting values :

4.7.1. Interface :

	- (void)readSettingValues;
Definition	Read BPM setting values

4.7.2. Delegate :

	- (void)WBOBLEManagerResponseReadSettingValue:(MicroLifeSettingValues *)settingValues;
Parameter	settingValues : BPM setting values

4.8. Write BPM setting values :

4.8.1. Interface :

	- (void)writeSettingValues:(MicroLifeSettingValues *)settingValues;
Definition	Write BPM setting values to device
Parameter	settingValues : BPM setting values

4.8.2. Delegate :

	- (void)WBOBLEManagerResponseWriteSettingValues:(BOOL)isSuccess;
Parameter	isSuccess : True or False

4.9. Read device ID and info :

4.9.1. Interface :

	- (void)readDeviceIDAndInfo;
Definition	Read device ID and info

4.9.2. Delegate :

	- (void)WBOBLEManagerResponseReadDeviceIDAndInfo:(MicroLifeDeviceInfo *)deviceInfo;
Parameter	deviceInfo : Device ID and info

4.10. Read device time :

4.10.1. Interface :

	- (void)readDeviceTime;
Definition	Read device time

4.10.2. Delegate :

	- (void)WBOBLEManagerResponseReadDeviceTime:(MicroLifeDeviceInfo *)deviceInfo;
Parameter	deviceInfo : Device Time

4.11. Write device Time

4.11.1. Interface :

	- (void)writeDeviceTime;
Definition	Write device Time to device

4.11.2. Delegate :

	- (void)WBOBLEManagerResponseWriteDeviceTime:(BOOL)isSuccess;
Parameter	isSuccess : True or False

4.12. Read BPM function setting value :

4.12.1. Interface :

	- (void)readFunctionSettingValue;
Definition	Read BPM function setting value

4.12.2. Delegate :

	- (void)WBOBLEManagerResponseReadFunctionSettingValue:(MicroLifeFunctionSettingValues *)functionSettingValues;
Parameter	functionSettingValues : BPM function setting value

4.13. Read BT module name :

4.13.1. Interface :

	- (void)readBTModuleName;
Definition	Read BT module name of device

4.13.2. Delegate :

	- (void)WBOBLEManagerResponseReadBTModuleName:(NSString *)BTModuleName;
Parameter	BTModuleName : BT Module Name

4.14 Start/Stop remote measurement :

4.14.1 Interface :

	- (void)startRemoteMeasurementWhithCBPFunction:(Dformat)d format;
Definition	Start remote measurement
Parameter	dformat : Data format which BPM sent out. It is same as what APP requested. NoCBPRaw : No CBP raw data LowCBPRaw : low resolution CBP data (sampling rate =16Hz) FullCBPRaw : full CBP raw data (sampling rate=256Hz)

	- (void)stopRemoteMeasurement;
Definition	Stop remote measurement

4.14.2 Delegate :

	- (void)WBOBLEManagerResponseStartRemoteMeasurement:(Dformat)D_format;
Definition	Start remote measurement
Parameter	dformat : Data format which BPM sent out. It is same as what APP requested. NoCBPRaw : No CBP raw data LowCBPRaw : low resolution CBP data (sampling rate =16Hz) FullCBPRaw : full CBP raw data (sampling rate=256Hz)

	- (void)WBOBLEManagerResponseremoteMeasurementStatus Every5secondsWithSTATUS:(STATUS)status MeasurementNumber:(int)measurementNumber TotalMeasurementNumber:(int)totalMeasurementNumber Countdown:(int)countdown TotalMeasuretime:(int)totalMeasuretime;
--	--


Definition	Send remote measurement status every 5 seconds
Parameter	<p>status :</p> <p>MeasurementWait : device wait countdown for next measurement.</p> <p>MeasurementStart : device is start BP measurement,</p> <p>MeasurementStop : manual press I/O to stop measurement</p> <p>measurementNumber : Send current measurement number in auto mode.</p> <p>totalMeasurementNumber : Send total measurement number in auto mode.</p> <p>countdown : Send current countdown time in auto mode.</p> <p>totalMeasuretime : Send total measurement time (seconds) in auto mode. Total measurement are count between 1st measurement to last measurement. (exclude rest time)</p>

	<p>-</p> <p>(void)WBOBLEManagerResponseMeasurementResultsForEachMeasurement:(MicroLifeCurrentAndMData *)data</p> <p>HistoryMeasurementNumber:(int)historyMeasuremeNumber</p> <p>CurrentMeasurementTimes:(int)currentMeasurementTimes</p> <p>AverageCalculationWhenMeasurement:(BOOL)isAverage;</p>
Definition	Send measurement results for each measurement
Parameter	<p>dRecord : CurrentAndMData</p> <p>historyMeasuremeNumber : The history measurement times store in memory.</p> <p>currentMeasurementTimes : Send current measurement times to APP.</p> <p>isAverage : Send Average calculation when measurement to APP.</p>


Chapter 5 User Interface & Functionality of WatchBP Office

5.1. Getting Started :

Start the app and then select the button “WatchBP Office II SDK” / “E” to communicate with the designate device.



5.2. Operating Interface and Sequence :



- 5.2.1. Region A : The log window is used to display information about communication handshake between App and device.
- 5.2.2. Region B : The display is for current device with MAC address.
- 5.2.3. Region C : This part is for all functionalities of WatchBP Office device.

C1 : Read CBP data by different mode such as Low/ Full. Before running this, the function “Read all history data” shall be performed firstly.

C2 : Read or Write User ID.

C3 : Read/ Write BPM parameters/ setting values :

1. Highest inflation pressure of AUS mode : Valid parameter: 0(not setting), 160, 180, 200, 220, 240
2. Highest inflation pressure of Auto mode : Valid parameter: 0(not setting), 160, 180, 200, 220,

240. The device inflates the cuff using fuzzy logic to proper cuff pressure; the Highest Inflation Pressure is considered as a safeguard pressure.
3. SW_AUTO_hide : Set Show readings during rest time in auto mode. true:hide/false:show.
 4. SW_SEL_silent : Beeper
 5. SW_AUS_Hide : Set Show cuff pressure during deflation in AUS mode. true:hide/false:show.
 6. SW_AVG_no_include_first : Set Average is include first memory data. true:is/false:is not.
 7. SW_CBP : Set CBP measurement
true:enabled/false:disabled.
 8. SW_AFib : Set AFib measurement
true:enabled/false:disabled.
 9. SW_AMPM : Set 12/24-hour clock true:12-hour/false:24-hour.
 10. SW_Kpa : Set Pressure unit:
true:Kpa/false:mmHg.
 11. RestTime : Rest time of auto mode. Start countdown base on rest time before 1st measurement in auto mode.
 12. IntervalTime : Interval time of auto mode. Start countdown base on interval time before 2nd~6th measurement in auto mode.
 13. AutoMeasureNumber: It's number of measurements in auto mode.

Note: Below are illegal settings

(1) AutoMeasureNumber=1 or 2 and


SW_AVG_no_include_first=true

(2) AutoMeasureNumber=1 and SW_AFib =true

(3) SW_CBP=true, AutoMeasureNumber *

IntervalTime >= 1800

For instance, the BPM setting values is included 13 parameters such as “160, 200, True, False, False, False, True, True, True, False, 15, 15, 3”



5.2.4. Refer to WBOViewController from the demo code (sample code) to get more detailed.

5.2.5. Operation Sequence :

5.2.5.1. The scanning (discovery) is automatically run to discover devices in the vicinity.


5.2.5.2. If a device is bonded, it will be connected accordingly. If not, the “bindingDevice” can be used to run bonding process.

5.2.5.3. Once the device is connected, select each function for communication with device.

Chapter 6 Demo App


Each “XXXBLEManager” has the same operation sequences and features. The following demonstration is mainly conducted by using the model WBPBLEManager. Refer to the sample code to see the complete application for dedicated devices.

6.1. Scanning / Discovery :


 <p>The screenshot shows the WBPBLEManager app interface. At the top, there's a status bar with the time 9:33 and signal indicators. Below the status bar is a navigation bar with a back arrow and the title 'WBPBLEManager'. The main content area displays a log of events with timestamps and messages. The log shows the app starting, the binding device being set to null, the Bluetooth state being updated to 'MicroLifeBLEStatePoweredOn', and a discovered Bluetooth device with MAC address <e8a7a342 99e3> and name 'WatchBP Home'. Below the log, there's a list of actions: getBindingDevice, bindingDevice, unBindingDevice, reScan, stopScan, and disconnectDevice. Each action has a corresponding description. At the bottom, there's a footer with the text 'DeviceSDK Version: 1.3'.</p>	<ol style="list-style-type: none"> 1. First of all, the WBPBLEManager starts Scanning / Discovery procedure. The detailed information will be displayed in Region A. 2. The WBPBLEManager bpmBindingDevice: * : (null) indicates that the remote device wasn't bonded yet. 3. The WBPBLEManagerCellPhoneBluetoothDidUpdateState:MicroLifeBLEStatePoweredOn indicates that the power status of Bluetooth from cell phone is On. 4. The WBPBLEManagerDidDiscoverBluetoothDeviceMacAddress:
---	--

	<p><e8a7a342 99e3> Name:WatchBP Home indicates that the device with MAC address <e8a7a342 99e3> is discovered nearby.</p> <p>5. In Region B, the “Did Discover Bluetooth Device : <e8a7a342 99e3>” stands for that the discovered/ selected device can be used in next steps such as Connection or Bonding.</p> <p>6. In Region C, click the function “bindingDevice” to run bonding process.</p>
--	---

6.2. Connection :

 <p>The screenshot shows the WBPBLEManager app interface. At the top, there's a status bar with the time 9:33 and battery level. Below it, a navigation bar with a back arrow and the title 'WBPBLEManager'. The main content area displays a log of events: <ul style="list-style-type: none"> te:MicroLifeBLEStatePoweredOn 2019-05-31 09:33:20:956 WBPBLEManagerDidDiscoverBluetoothDeviceMac Address:<e8a7a342 99e3> Name:WatchBP Home 2019-05-31 09:33:32:724 WBPBLEManager bindingDevice:<e8a7a342 99e3> 2019-05-31 09:33:33:451 WBPBLEManagerDidConnectDevice:{ RSSI = "-59"; UUID = "E24391D9-FE56-33C6-4871-DCB81A306EC7"; macAddress = <e8a7a342 99e3>; name = "WatchBP Home"; } Below the log, there are several buttons with corresponding actions: <ul style="list-style-type: none"> getBindingDevice (Read usual mode history data) bindingDevice (Read diagnostic mode history data) unBindingDevice (Clear selected mode history data) reScan (Clear current mode history data) stopScan (Disconnect the Bluetooth) disconnectDevice (Write device Time) (Write a new user ID) At the bottom, there's a footer with 'DeviceSDK Version: 1.3'. </p>	<p>1. The information about the device bonding will be shown “WBPBLEManager bindingDevice:<e8a7a342 99e3>” in Region A.</p> <p>2. The following indicates that the device is connected. WBPBLEManagerDidConnect Device:{</p> <pre> RSSI = "-68"; UUID = "E24391D9-FE56-33C6-4871-DCB81A306EC7"; macAddress = <e8a7a342 99e3>; name = "WatchBP Home"; </pre> <p>}</p>
--	--

6.3. Bonding / Pairing :



The screenshot shows the WBPBLEManager app interface. At the top, there is a status bar with the time 10:21 and signal indicators. Below the status bar, there is a navigation bar with a back arrow and the title WBPBLEManager. The main content area displays a log of events, including a Bluetooth pairing request from "WatchBP Home". A modal dialog box titled "Bluetooth Pairing Request" is overlaid on the log, asking the user to confirm the pairing with the code "517583". The dialog has two buttons: "Cancel" and "Pair". Below the dialog, there is a list of actions with their corresponding descriptions: "getBindingDevice" (Read usual mode history data), "bindingDevice" (Read diagnostic mode history data), "unBindingDevice" (Clear selected mode history data), "reScan" (Clear current mode history data), "stopScan" (Disconnect the Bluetooth), "disconnectDevice" (Write device Time), and "Write a new user ID". At the bottom, there is a footer with the text "DeviceSDK Version: 1.3".


1. There is a message to confirm the bonding between device and cellphone if they haven't bonded yet.

2. Click "Pair" to run the bonding process.

6.4. Testing Command : Write a new user ID to BPM


<div data-bbox="499 331 906 365"> < Back WBPBLEManager </div> <div data-bbox="515 398 898 745"> <pre> 6854B8F8C320906854B9090420906854B91145 20906854B9198620906854B921C720906854B9 2A0821906854B9324920916955B93A8A213D *** 2019-06-11 10:24:46:405 WBPBLEManager Write a new user ID: 835702427AY *** 2019-06-11 10:24:46:603 4DFF18060038333537303234323741590000000 000000000000DA *** 2019-06-11 10:24:46:931 WBPBLEManagerResponseWriteUserID:YES *** 2019-06-11 10:24:47:110 4D5103068128 *** </pre> </div> <div data-bbox="515 779 898 801"> Did Discover Bluetooth Device: <e8a7a342 99e3> </div> <div data-bbox="515 813 898 1081"> <div> getBindingDevice Read usual mode history data </div> <div> bindingDevice Read diagnostic mode history data </div> <div> unBindingDevice Clear selected mode history data </div> <div> reScan Clear current mode history data </div> <div> stopScan Disconnect the Bluetooth </div> <div> disconnectDevice Write device Time </div> <div> Write a new user ID </div> </div> <div data-bbox="675 1115 738 1137"> 11 11 11 11 </div>	<ol style="list-style-type: none"> 1. The command “Write a new user ID” is to write a new user ID to device. 2. The following “WBPBLEManager Write a new user ID : 835702427AY” indicates that the writing value is “835702427AY” made up of ASCII code. 3. The response is as below: WBPBLEManagerResponse WriteUserID:YES. It means that the process is successful. 4. It can be checked by the command “Read user ID and version data”.
---	--

6.5. Testing Command : Read user ID and version data from BPM


 <p>2019-06-12 11:23:46:000 WBPBLEManagerResponseReadUserAndVersionData: ta: BPMUserID = 500272727ON; EBodyuserID = "<null>"; UUserID = "<null>"; age = 0; gender = 0; height = 0; resistance = 0; roleType = 0; userNo = 1; weight = 0; } VersionData: FWName = rD7; currentMode = 0; day = 18; deviceBatteryVoltage = 0; errHappendTimes = "<null>";</p> <p>Did Discover Bluetooth Device: <e8a7a342 99e3></p> <p>getBindingDevice Read nocturnal mode setting bindingDevice Change nocturnal mode setting unBindingDevice Erase all of device measure & error times reScan Read device ID and info stopScan Read device Time disconnectDevice Read user ID and version data disconnectDevice Read Nocturnal mode history data</p> <p>9:41</p> <p>WBPBLEManager</p> <p>} VersionData: FWName = rD7; currentMode = 0; day = 18; deviceBatteryVoltage = 0; errHappendTimes = "<null>"; hour = 0; isTimeReady = 0; macAddress = "<null>"; maxMemory = 250; maxUser = 1; measurementTimes = 0; minute = 0; month = 4; openBPtype = 0; openNoUsualModeAFib = 0; openNocturnal = 0; openNocturnalMode = 0; optionAFib = 0;</p> <p>Did Discover Bluetooth Device: <e8a7a342 99e3></p> <p>getBindingDevice Read nocturnal mode setting bindingDevice Change nocturnal mode setting unBindingDevice Erase all of device measure & error times reScan Read device ID and info stopScan Read device Time disconnectDevice Read user ID and version data disconnectDevice Read Nocturnal mode history data</p> <p>Label</p>	<p>1. The command "Read user ID and version data" is to receive the information of user ID and device version.</p> <p>2. The response is as below: WBPBLEManagerResponseReadUserAndVersionData: BPMUserID = 500272727ON; EBodyuserID = "<null>"; UUserID = "<null>"; age = 0; gender = 0; height = 0; resistance = 0; roleType = 0; userNo = 1; weight = 0; } VersionData: FWName = rD7; currentMode = 0; day = 18; deviceBatteryVoltage = 0; errHappendTimes = "<null>"; hour = 0; isTimeReady = 0; macAddress = "<null>"; maxMemory = 250;</p>
---	--

	<pre>maxUser = 1; measurementTimes = 0; minute = 0; month = 4; openBPtype = 0; openNoUsualModeAFib = 0; openNocturnal = 0; openNocturnalMode = 0; optionAfib = 0; optionAmbientT = 1; optionDeviceID = 0; optionDiagnosticModeAFib = 0; optionIHB = 1; optionMAM = 1; optionTubeless = 0; protocolVersion = 0; second = 0; year = 2019; }</pre> <p>3. The “BPMuserID” stands for user ID and the value is “500272727ON”.</p> <p>4. The “VersionData” stands for device version included name of firmware “FWName” and its building date “day”, “month” and “year”.</p>
--	--


6.6. Disconnection

 <p>The screenshot shows the WBPBLEManager app interface. At the top, the status bar displays the time 9:48 and signal icons. Below the title bar, there is a log of Bluetooth operations. The log includes a long MAC address, a timestamp 2019-05-31 09:46:43:993, and the message 'WBPBLEManager Disconnect the Bluetooth'. Below this, another timestamp 2019-05-31 09:46:44:194 and the MAC address 4DFF020452 are shown. Further down, a timestamp 2019-05-31 09:46:45:107 and the message 'WBPBLEManagerDidFailToConnectDevice' are displayed. Below the log, there is a message 'Did Discover Bluetooth Device: <e8a7a342 99e3>'. At the bottom, there is a list of actions with their corresponding descriptions: 'getBindingDevice' (Read usual mode history data), 'bindingDevice' (Read diagnostic mode history data), 'unBindingDevice' (Clear selected mode history data), 'reScan' (Clear current mode history data), 'stopScan' (Disconnect the Bluetooth), and 'disconnectDevice' (Write device Time). At the very bottom, the text 'DeviceSDK Version: 1.3' is displayed.</p>	<p>1. In order to disconnect the device completely, either “disconnectDevice” (Region C) or “Disconnect the Bluetooth” (Region D) is available.</p> <p>2. The response is the following “WBPBLEManagerDidFailToConnectDevice” once the disconnection command is executed.</p>
---	---

6.7. Testing Command : Read BPM setting values (WBOBLEManager)

 <p>The screenshot shows the WBOBLEManager app interface. At the top, the title is 'WBOBLEManager'. Below it, there's a log area showing the command 'WBOBLEManager Read BPM setting values' and its response. The response is a JSON object containing various BPM settings. Below the log, there are several interactive elements: a 'Read all history data' button, a 'Clear all history' button, a 'Read CBP data' button, a 'Read user ID and version' button, a 'Write a new user ID' button, a 'Read BPM setting values' button, a 'Write BPM setting values' button, a 'Read device Time' button, a 'Write device Time' button, a 'Read BPM function setting value' button, a 'Read BT module name' button, and a 'Read device ID and info from BPM' button. At the bottom, there's a 'Label' field with a placeholder text 'Label'.</p>	<p>1. The command “Read BPM setting values” is to read the setting values of BPM.</p> <p>2. The response is as below: WBOBLEManagerResponseReadSettingValue:{</p> <pre> ABPMEnd = 0; "ABPMEnt_first" = 0; "ABPMEnt_second" = 0; ABPMStart = 0; "AUS_HI_infPressure" = 160; AutoMeasureNumber = 3; "CBPInt_first" = 0; "CBPInt_second" = 0; "CBP_zone1_meas_off" = 0; "CBP_zone2_meas_off" = 0; "HI_infPressure" = 200; IntervalTime = 15; RestTime = 15; "SW_AFib" = 1; "SW_AMPM" = 1; "SW_AUS_Hide" = 0; "SW_AUTO_hide" = 1; "SW_AVG_no_include_first" = 0; "SW_CBP" = 1; "SW_Kpa" = 0; "SW_SEL_silent" = 0; "SW_checkhide" = 0; } </pre>
--	--

6.8. Testing Command : Start remote measurement / Stop remote measurement (WBOBLEManager)

	<p>1. The command “Start remote measurement” / “Stop remote measurement” is to run the process of remote control.</p> <p>2. During measuring process, the relative status can be found by the following response:</p> <pre>WBOBLEManagerResponseRemoteMeasurementStatusEvery5secondsWithSTATUS:1 MeasurementNumber:1 TotalMeasurementNumber:3 Countdown:15 TotalMeasuretime:0</pre> <p>3. Once the measurement is completed, the result is the following:</p> <pre>WBOBLEManagerResponseMeasurementResultsForEachMeasurement:{ AA = 0; ABPM = 0; AFib = 0; AM = 0;</pre>
--	--

	<pre>CPP = 60; CSBP = 111; IHB = 0; LB = 0; MAM = 0; MAP = 58; MCBP = 786; PM = 0; PVR = 202; SM = 0; arr = 0; condition = 1; cuffokr = 0; day = 26; deviceMode = 66; dia = 52; diagnostic = 0; errorCode = 0; errorMessage = "<null>"; hour = 16; hr = 76; isFor3G = 0; minute = 27; month = 11; resultCode = 0; sys = 103; usual = 0; year = 2019; } HistoryMeasurementNumber: 0 CurrentMeasurementTimes:3</pre>
--	--

	<p>AverageCalculationWhenMeasurement: YES</p> <p>◦</p> <p>3.2.WBOBLEManagerResponseReadCBPData:{</p> <p> CBPdatas = (</p> <p>);</p> <p> CalCBPs = (</p> <p> 5469,</p> <p> 5461,</p> <p> 5446,</p> <p> </p> <p> 5492</p> <p>);</p> <p> dformat = 0;</p> <p>}</p>
--	--